

FLJ16126-hypothetical LOC645010 Gene

Reference: PTXBC130506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC645010
GENE ID	645010
ORIGIN SPECIES	Human
ACCESSIONS	BC130506
VECTOR	pDONR223
NCBI SYMBOL	FLJ16126
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattgcacattctgtccaaccagagtcctgctcctaagcacagcatcttctctgttaattatcccttttaacagcaacaaatgacctggccggcacacagacacagatgagggccacacctgtaattctgggtctgagggggaccctgggtggttctgggtgctg
cgggttcttcagtaaaagaacttcaggaggttctaatcagctgtttcatcaagattctcataacctctctgctgccgtttcaagtgggtttgttcccgagcacattattaggaaagcggctataaactgttacctgcctccggctccgctccacaacacccgccaagcccacactgtgca
aaacagtag