

LELP1-late cornified envelope-like proline-rich 1 Gene

Reference: PTXBC130507

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	late cornified envelope-like proline-rich protein 1; novel small proline rich protein; late cornified envelope like proline rich 1		

For research use only.

GENE DESCRIPTION	late cornified envelope-like proline-rich 1
GENE ID	149018
ORIGIN SPECIES	Human
ACCESSIONS	BC130507
VECTOR	pDONR223
NCBI SYMBOL	LELP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagtgatgataaaagtaaatcaaatgacccaagactgagcccaagaactgcgatccaagtgtgaacaaaagtgtgagtcgaatgccagccagctgtttaagaagctgctgcaacgctgttcgaaaagtgccatgggaaaagtgtccagcaccaccaagtgctgcct
gccccctgcagtccttcctcctgcccctccccagccctgcaccaagccctgtcctcctaaatgcccttcctcctgcccacatgcttgcacacctcctgcctccccagagtga