



Reference: PTXBC130529

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HAKAI; RNF188; E3 ubiquitin-protein ligase Hakai; Cas-Br-M (murine) ecotropic retroviral transforming sequence-like 1; Cbl proto-oncogene like 1, E3 ubiquitin protein ligase; Cbl proto-oncogene, E3 ubiquitin protein ligase-like 1; E-cadherin binding prote		

For research use only.

GENE DESCRIPTION	Cas-Br-M (murine) ecotropic retroviral transforming sequence-like 1
GENE ID	79872
ORIGIN SPECIES	Human
ACCESSIONS	BC130529
VECTOR	pDONR223
NCBI SYMBOL	CBLL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1