

TR2IT1-thioredoxin reductase 2 intronic transcript 1 Gene

Reference: PTXBC130546

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TR2IT1; TXNRD3IT1; TXNRD3NT1; protein TXNRD3NB; TXNRD3 neighbor gene protein; thioredoxin reductase 2 intronic transcript 1; thioredoxin reductase 3 intronic transcript 1; thioredoxin reductase 3 neighbor gene protein; thioredoxin reductase 3 new transcri		

For research use only.

GENE DESCRIPTION	thioredoxin reductase 2 intronic transcript 1
GENE ID	645840
ORIGIN SPECIES	Human
ACCESSIONS	BC130546
VECTOR	pDONR223
NCBI SYMBOL	TR2IT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggacacgagcagagagagacttggcagccggagctgaaagctgagcagcagatgccactggagcctgtgagggcacgattgagtggtggcctgctgctctccacaccactgcagaggcgtcctctctggaacatggagacaagggtttggacaggggtttccca
gtccctgggaagagattaagagactactgaagatccagggcactgcaagccagatctgtacctccaccagagagaaggcaaatgtcttctggagagcctggcagccagaggggaaagggtcaagaaacctatcctggccctgggaaagtggagggttaaggcagagccagcca
tgagaagagatgatgttggccagcgaatgtatcagcggatga