

MUSTN1-musculoskeletal, embryonic nuclear protein 1 Gene

Reference: PTXBC130550

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MUSTANG; musculoskeletal embryonic nuclear protein 1; musculoskeletal temporally activated novel; musculoskeletal, embryonic nuclear protein 1		

For research use only.

GENE DESCRIPTION	musculoskeletal, embryonic nuclear protein 1
GENE ID	389125
ORIGIN SPECIES	Human
ACCESSIONS	BC130550
VECTOR	pDONR223
NCBI SYMBOL	MUSTN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaggctgtgtcaggaagccccctatcaagaagaagcggccccctgtgaaggaggaggaacctgaagggggcccgaggaaacctgaccaagaaccaggaaatcaagtccaagacctaccaggtcatgcgagagtgtgagcaagctggctcggccggccccgtcggtgttca
ggcgacccgcacaggtaccgagactgtctttgagaagcccaagccggacccaccaagagtgtcttcggctga