



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CCL25-chemokine (C-C motif) ligand 25 Gene

Reference: PTXBC130561

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

Ckb15; SCYA25; C-C motif chemokine 25; Ck beta-15; TECKvar; chemokine (C-C motif) ligand 25; chemokine TECK; small inducible cytokine subfamily A (Cys-Cys), member 25; small-inducible cytokine A25; thymus expressed chemokine; C-C motif chemokine ligand 25

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 25
GENE ID	6370
ORIGIN SPECIES	Human
ACCESSIONS	BC130561
VECTOR	pDONR223
NCBI SYMBOL	CCL25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacctctggctctggcctgcctggctggccgctctctcggagcctggcccccctgctccacacccaaggctgcttctgaggactgctcctggcctaccactaccctatgggtggcgtgctccggcgccgctgacttaccgagatccagagagtgagcgggagctgcaatctg
cctctgctgatatttaccctccccagagacagaggaagctgctggcagaccccaaaagcagggagctgcagagagccatgaagctctctgcatgctcgaataaaggctttgcaaaagctccacacaacacgcagacctccaagacccctcatgctgtaaaagattgagttctgga

More info: Please contact us for any question: contact@proteogenix.fr