

PLAC8L1-PLAC8-like 1 Gene

Reference: PTXBC130586

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLAC8-like protein 1; PLAC8 like 1		
For research use only.			

GENE DESCRIPTION	PLAC8-like 1
GENE ID	153770
ORIGIN SPECIES	Human
ACCESSIONS	BC130586
VECTOR	pDONR223
NCBI SYMBOL	PLAC8L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattggttggaaagtaacttctcaggtgtcctgaggatcttcttactcaacataatactacacctctgttgcacacatgtcatctgaagatgaacattttattccaacttgagagccatgtaccagccagcgtgtggtgaagcagcctgttcggggagccagtggcaggacgacaatcac
agcaattgtccagactggcggggctggagcaccggtctctcaggtgtcgcagagatggagaattgttctgtgtctattctgtcctatgtgtcttgagtgatcgccaggcattatggagagtgtctttgtggccgtgtgttacctgggtccaccttgcactgagaattggcaccaggg
agagacataaaatacaggggcacactgtgtgaagactggctggcgggtgcactgctgttgggctttttccatctgccaggtggccgggaactcaagatgaggacctccaagtctatgaaatctgtgcagtccaatgactaaggacacctggttga