

POM121C-POM121 membrane glycoprotein C Gene

Reference: PTXBC130587

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POM121-2; nuclear envelope pore membrane protein POM 121C; POM121 membrane glycoprotein C; nuclear pore membrane protein 121-2; pore membrane protein of 121 kDa C; POM121 transmembrane nucleoporin C		

For research use only.

GENE DESCRIPTION	POM121 membrane glycoprotein C
GENE ID	100101267
ORIGIN SPECIES	Human
ACCESSIONS	BC130587
VECTOR	pDONR223
NCBI SYMBOL	POM121C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacacattgcaggggccagtgatcattcaagatgtggctgtggttcaccaggaggagtgccgcaactggaccctgatgagaagataacacacgggatgtgatgtggagaactacagccatctagttccttgcttatgaggtggcaacatctgtacttcggagattctgaagc
cgagcaacttgcccaagtctcttctttcccattaa