

C10orf114-chromosome 10 open reading frame 114 Gene

Reference: PTXBC130594

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf114; bA418C1.3; protein CASC10; cancer susceptibility candidate gene 10 protein; cancer susceptibility candidate 10		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 114
GENE ID	399726
ORIGIN SPECIES	Human
ACCESSIONS	BC130594
VECTOR	pDONR223
NCBI SYMBOL	C10orf114
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaatcggggagccctcagggtggcgacggggagaggcagcggggtggcgctgcagaggcgtccccacccctccggggaccggggcccgaggctgccaggcccgtgcacgcccggcgccgggagaaacacggggcagcagccccctgcaggtgccgg
gcgcctccgcagccgtgcgcggaccgactcttaagggtggcaccacgggttcccagcccgcgcaaccgcagccctggctccatccgcccacttctccctgtagcggcgccctgaccgtccggaaacggccggagtgctgatgcctgctgtatctgctgatccctcacc
ctacctcatccaggactcttccgcgatgcgcagtagcgaatttcagagcgagattga