

C21orf34-chromosome 21 open reading frame 34 Gene

Reference: PTXBC130602

Product type: DNA Storage Condition: -20°C
Size/Quantity: 2 ug Delivery Condition: +4°C, Blue Ice
Synonyms: C21orf34; C21orf35; FLJ38295; long intergenic non-protein coding RNA 478
For research use only.

GENE DESCRIPTION chromosome 21 open reading frame 34
GENE ID 388815
ORIGIN SPECIES Human
ACCESSIONS BC130602
VECTOR pDONR223
NCBI SYMBOL C21orf34
LENGTH ProteoGenix bp
SEQUENCE PRIMERS Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT
Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaattatattagaacttcagcaaatgacagctgcataatagagaagccgactcagattgacagccaataagccagctgtggaagggcctctgggtgacagctctggcttttgggtgttgacaacctagggcaaaactgaaatggacagggcagaaaatagcacac
ttttactctacgaactgttcttttcattgatctcaggaattttactggataattatccaataaattgcaagcattctccaaatggagctctttctgagatgaagagaattctcaatgcaagattga