



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LCE2D-late cornified envelope 2D Gene

Reference: PTXBC130603

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

LEP12; SPRL1A; late cornified envelope protein 2D; late envelope protein 12; small proline rich-like (epidermal differentiation complex) 1A; small proline-rich-like epidermal differentiation complex protein 1A; late cornified envelope 2D

For research use only.

GENE DESCRIPTION	late cornified envelope 2D
GENE ID	353141
ORIGIN SPECIES	Human
ACCESSIONS	BC130603
VECTOR	pDONR223
NCBI SYMBOL	LCE2D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]