

SCGB1D4-secretoglobin, family 1D, member 4 Gene

Reference: PTXBC130639

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IIS; secretoglobin family 1D member 4; IFN-gamma inducible SCGB (IIS); IFN-gamma-inducible secretoglobin		

For research use only.

GENE DESCRIPTION	secretoglobin, family 1D, member 4
GENE ID	404552
ORIGIN SPECIES	Human
ACCESSIONS	BC130639
VECTOR	pDONR223
NCBI SYMBOL	SCGB1D4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctgcagtggtctcctgatggctcgtggcccttgcgtctaccaggcccatgctctgtctgccagctgttgccttgagatcacagctcttattcttaagtgcgctgcggttaacctcaagttgccaaactaatccacctcagaagctcttcagccaagttggaagtgaagcactgcaccgatcagatatcttttaagaacgactctcattgaaaaagtcctggtggaaatag