

BLOC1S1-biogenesis of lysosomal organelles complex-1, subunit 1 Gene

Reference: PTXBC130640

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BLOS1; BORCS1; GCN5L1; MICoA; RT14; biogenesis of lysosome-related organelles complex 1 subunit 1; BLOC subunit 1; BLOC-1 subunit 1; GCN5 (general control of amino-acid synthesis, yeast, homolog)-like 1; GCN5 general control of amino-acid synthesis 5-like		

For research use only.

GENE DESCRIPTION	biogenesis of lysosomal organelles complex-1, subunit 1
GENE ID	2647
ORIGIN SPECIES	Human
ACCESSIONS	BC130640
VECTOR	pDONR223
NCBI SYMBOL	BLOC1S1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtccgcctcctaaagaacaccaggccaagcagaatgaacgcaaggagctgcaggaaaaggaggcgagaggctatcactgcagcgacctgcctgacagaagcttgggtgatcacctcaatgtgggtgtggccaggcctacatgaaccagagaaaagctggaccatga
gggtgaagaccctacaggtccaggctgcccaatttgcgaagcagacaggccagtggatcgggaatggtggagaactcaaccaggcactcaaggaaattggggatgtggagaactgggctcggagcatcgagctggacatgcgcaccattgccactgcaatgtctacaaagg
cagctgcagctgtgcccttcctag