



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM211-transmembrane protein 211 Gene

Reference: PTXBC130647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bA9F11.1; transmembrane protein 211		
For research use only.			

GENE DESCRIPTION	transmembrane protein 211
GENE ID	255349
ORIGIN SPECIES	Human
ACCESSIONS	BC130647
VECTOR	pDONR223
NCBI SYMBOL	TMEM211
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctctcggagctggtctctgttggccctcaatgcaatttctctgtcttgggctgtggcccccaaggctgtgcccaaggagaagcagttgtccaatgccaggggtgcaggcagttggcagctactgccatgattgtgggtctgctgattttccaatcggcctgtcctcccaatca
caaggaaagtgtcgaagctctctccatgtattatgtgggaagtgcggcgtgttgggtttacatgctctatctcaatcgactctgtgccagcctctgccccatcatcagctgtccccacacaacaaggctcgaaggagagaccatcatcttccagtgccaccgagagaatcatct
ttgtcccaacattgaacaattaa