

KU-MEL-3-KU-MEL-3 Gene

Reference: PTXBC130651

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	KU-MEL-3
GENE ID	497048
ORIGIN SPECIES	Human
ACCESSIONS	BC130651
VECTOR	pDONR223
NCBI SYMBOL	KU-MEL-3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgcgaactcctgcgccgacatcctctccagctcctacctgcaggcctcctgacctcaggctcctaccctccgctggcctctgtctacactcatctttgtgcctgcttccaatcatgccaggcctggggctcctgtcttaactcaccatga