

PF4V1-platelet factor 4 variant 1 Gene

Reference: PTXBC130653

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXCL4L1; CXCL4V1; PF4-ALT; PF4A; SCYB4V1; platelet factor 4 variant; C-X-C motif chemokine 4; PF4alt; PF4var1; platelet factor 4, variant 1 (PF4-like); platelet factor 4 variant 1		

For research use only.

GENE DESCRIPTION	platelet factor 4 variant 1
GENE ID	5197
ORIGIN SPECIES	Human
ACCESSIONS	BC130653
VECTOR	pDONR223
NCBI SYMBOL	PF4V1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctccgagccaggtccgctcaccgcgccaccgccagagatgctgttctggcgttgctcctgccagttgtggtcgccttcgccagagctgaagctgaagaagatggggacctgcagtgccctgtgtgaagaccacctcccaggtccgtcccaggcacatcaccagccfaggagtgatcaagccggaccccactgcccactgcccactcatagccacgtgaagaatggaggaaaatttgctggatctgcaagccctgctgtacaagaaaatcataaggaacatttgagaggttag