

LOC642947-hypothetical protein LOC642947 Gene

Reference: PTXBC130654

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC642947
GENE ID	642947
ORIGIN SPECIES	Human
ACCESSIONS	BC130654
VECTOR	pDONR223
NCBI SYMBOL	LOC642947
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggccaggcagggactgtggaacaggcaaccaggatgggggttggggggggcgctcaaatctgacaagccctgtgccatgggcaagatgtctgtctgtctaggtccaatagtcagcaaaagaccaaagctacctagagggtgtgtgtgagctttgggagatggacatacctgg
ccatgtctccactgcagccattccacaccaggcccttctgggtccgtgcagactggagtcctgtttctgccaacttccaagcagctgtctttgcaacctcaaatgtctgtagggtcatggagcctcttcagctaggattctagaggccatgacaagaacaggctactccatgctgtttaac
tcatcctttcccaggagcatttgggtcagaatgagtcagtggtgctag