

C12orf59-chromosome 12 open reading frame 59 Gene

Reference: PTXBC131523

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf59; transmembrane protein 52B		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 59
GENE ID	120939
ORIGIN SPECIES	Human
ACCESSIONS	BC131523
VECTOR	pDONR223
NCBI SYMBOL	C12orf59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagtcgagttcatgtcgtggcgccctcagccctgtgtatttcacctgtttctgggacgagatgtgaggaaactgtgtaactcctgaacattgcctgaccacagactgggtacatctctggatatatggttctagtggttaattggcgcgctgttctcctgtgtggcctgacgtccc
tggcttccgctgtgtgtgtgagccgagcaaaatggggaagatggggccaccacccctgtgaagtgcctgtcattgttcgatcagcagcactctccagagcactatcacatctctgcagtcgggtgttggccctgcagctcggagatcctggctgtgtcactccacagct
ccctgggcccagctgcctcctcttggacacctcccagggtatgaagaagctcttcacatgagtcgttcacagtagccatgtgcgggcagaaagcacctgatctaccccgagctacctgaagaaagcagctgcctccaacagagaaggagtcgactcgaatagttgactcttgaactg
a