



Reference: PTXBC131524

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LSDP5; LSDA5; MLDP; OXPAT; perilipin-5; lipid storage droplet protein 5; perilipin 5		

For research use only.

GENE DESCRIPTION	lipid storage droplet protein 5
GENE ID	440503
ORIGIN SPECIES	Human
ACCESSIONS	BC131524
VECTOR	pDONR223
NCBI SYMBOL	LSDP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1