

TIAF1-TGFB1-induced anti-apoptotic factor 1 Gene

Reference: PTXBC131561

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAJN; SPR210; TGFB1-induced anti-apoptotic factor 1; 12 kDa TGF-beta-1-induced antiapoptotic factor; TGF-beta-1-induced antiapoptotic factor 1; molecule associated with Jak-3 N-terminal		

For research use only.

GENE DESCRIPTION	TGFB1-induced anti-apoptotic factor 1
GENE ID	9220
ORIGIN SPECIES	Human
ACCESSIONS	BC131561
VECTOR	pDONR223
NCBI SYMBOL	TIAF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcatctccatcttacccttcagagagcagtccttctctgtgcagctggagagcgtgtgagagagccgggtccaggttcttaagaatgaggtgaggaggggctctccggtgctgctgggctgggttgagcaagcctacgcagacaagtgtgtgtgtggaccatccgcacctccagc
ccccacccacccctcttctcagcgtgttatgtgcaatgacctatftaaggtaaacccattccaactacagcagttcagggctgatccaagcactgcctccctcctgctgtgtccaggtgtgtgtgacataaaactcaacttgagagggaaggcttgggttga