



Reference: PTXBC131580

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR147; NPFF1; NPFF1R1; OT7T022; neuropeptide FF receptor 1; G-protein coupled receptor 147; RFamide-related peptide receptor OT7T022; neuropeptide FF 1		

For research use only.

GENE DESCRIPTION	neuropeptide FF receptor 1
GENE ID	64106
ORIGIN SPECIES	Human
ACCESSIONS	BC131580
VECTOR	pDONR223
NCBI SYMBOL	NPFFR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]