

OR11A1-olfactory receptor, family 11, subfamily A, member 1 Gene

Reference: PTXBC131608

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	6M1-18; OR11A2; dJ994E9.6; hs6M1-18; olfactory receptor 11A1; olfactory receptor 11A2; olfactory receptor OR6-30; olfactory receptor, family 11, subfamily A, member 2; olfactory receptor family 11 subfamily A member 1		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 11, subfamily A, member 1
GENE ID	26531
ORIGIN SPECIES	Human
ACCESSIONS	BC131608
VECTOR	pDONR223
NCBI SYMBOL	OR11A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaattgtctccacaggaacgaaactattactgaattgtctccttggtctctatgacatccctgaactgcattctgtttttattgtattcactgctgtctatgctctcatcataggggaatatgctgattattgtagcagtggttagctccagaggtccacaaacctgtatatcttctgg
cgaatctgtccttctggaattttctacactccgcagtgatgccaaaatgctggagggtctctgcaagaagcaactatctctgtggctggttgcctgctccagttcttctctcggtctctagccacagctgaatgcttactgctgctgctcatggcatatgaccgtacctggcaatttgta
cccactccactaccactcctgatggggccagacggtacatggggctgggtgcacaacctggctctctggattgtggttagatggactgggttggccctggtgaccagctgaggttctgtggcccaaccacattgaccagttttactgtgactttatgcttttctgtggcctggctgctc
ggatcccagagtggtcaggtgacaactctcatctgtctgtgtctgctcactatcttttggactgattctgacatcttatccagaattgtggtggcagtgctgagagttcctgtggggcaagcaggagaagggtcttccacatgctctccacactagctgtgtgaccacattctatg
gaacgctcatgatctttatgtgaccctctgctgccatccagctcctccaaggtcttctcctgctctacactgtggtcaccctctctcaatcctgtgactataccatgaggaacaaggaggtgcatcaggcacttcggaagattctctgtatcaaaacactgaaacacttgattga