



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C2orf83-chromosome 2 open reading frame 83 Gene

Reference: PTXBC131618

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	folate transporter-like protein C2orf83; chromosome 2 open reading frame 83		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 83
GENE ID	56918
ORIGIN SPECIES	Human
ACCESSIONS	BC131618
VECTOR	pDONR223
NCBI SYMBOL	C2orf83
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaa gactat gccctg acgttt ggaaat taacct ccttcatt gctctg atgatt caaccat cgtgac gatgact gtgttg gacaac caaggact ggggtt cctctgttgacatt cagcaca aagcaagccctt ccacaa agccagccag atgtct gccact gctatcttg gacctt cccagcttcc aagcttcc aagctat taatg tccatct gaga aactg agcgcctt ggagacca gggatgag gacacatt ctagc tacaag catgctggcgtt gacttccatt cacttccctt cggggc cctctg cctatccact cagcagg ctctcagc cctcagc ggaactctctc cagcactg cctctaa acagacag atgtg agggag gacatg catgta gaaagg tccacag gaacac caagat ggtg cctctg cgtcag cactctag atgtg gcttctcagg ggaactag