



Reference: PTXBC131624

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BASS2; KRBOX2; LST003; pogo transposable element with KRAB domain; KRAB box domain containing 2		

For research use only.

GENE DESCRIPTION	pogo transposable element with KRAB domain
GENE ID	57645
ORIGIN SPECIES	Human
ACCESSIONS	BC131624
VECTOR	pDONR223
NCBI SYMBOL	POGK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]