



Reference: PTXBC131633

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD repeat-containing protein 27; WD repeat domain 27		
For research use only.			

GENE DESCRIPTION	WD repeat domain 27
GENE ID	253769
ORIGIN SPECIES	Human
ACCESSIONS	BC131633
VECTOR	pDONR223
NCBI SYMBOL	WDR27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr