

TMEM213-transmembrane protein 213 Gene

Reference: PTXBC131697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 213		
For research use only.			

GENE DESCRIPTION	transmembrane protein 213
GENE ID	155006
ORIGIN SPECIES	Human
ACCESSIONS	BC131697
VECTOR	pDONR223
NCBI SYMBOL	TMEM213
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgtttccccaacacccagcataagttcacctctggcagagttgtttccagctcctgcaggtgggagtcgactcactgcagcaggcactcggcacaactcgcaggaccggctcacctgcacggggcactcagcacagctccagcatgcagcgcctcccgcctgccaccg
ggccaccctgatctcagcctggcctttgctccctccactcggcctgtctcggcagcagcagcaacagctcaagcttgaccgtcaccaccagacccctgggacccctggagcagtgccctaacgtggacttctgcccacaagcagccgggtgctgccgacagagtgagcag
tacggctggatcgcggcagctgttgctggagcctctggttctcacctcatcctgtctgtgtggacaaactgatgaagctgactccagatgagccaaggacttgcagcgtga