



Reference: PTXBC131713

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanylate-binding protein 6; GBP-6; GTP-binding protein 6; guanine nucleotide-binding protein 6; guanylate binding protein family member 6		

For research use only.

GENE DESCRIPTION	guanylate binding protein family, member 6
GENE ID	163351
ORIGIN SPECIES	Human
ACCESSIONS	BC131713
VECTOR	pDONR223
NCBI SYMBOL	GBP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]