

C5orf51-chromosome 5 open reading frame 51 Gene

Reference: PTXBC131721

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0600 protein C5orf51; chromosome 5 open reading frame 51		
For research use only.			

GENE DESCRIPTION	chromosome 5 open reading frame 51
GENE ID	285636
ORIGIN SPECIES	Human
ACCESSIONS	BC131721
VECTOR	pDONR223
NCBI SYMBOL	C5orf51
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggccgagctctctgtgtggtgagacgagtggaagagctcgggagctcagggccacacagcaacttagcgaagctgccggtgaagatgatcacttttaattcgggcctctgcagcattagaaaaattgaaactcctgtgtggagaagagaagaatgttcaaatccatc
aaatctctagaactttacacacaggctattttgacatgacatattttgaggagaacaagctagtagatgaagattttcctgaagactctcttcacagaaagtaaaagagctgattttttcttcagaaccagaaatttttagtaaggaaaaataatgatccaaaacactgcaatttgcctggg
gatgagctactggaatgtctctccatttccgctgtggctgagccaaacttcatttctttaggcctctgcctccctag