

CAPS2-calcyphosine 2 Gene

Reference: PTXBC131724

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UG0636c06; calcyphosin-2; calcyphosphine 2; calcyphosine 2		
For research use only.			

GENE DESCRIPTION	calcyphosine 2
GENE ID	84698
ORIGIN SPECIES	Human
ACCESSIONS	BC131724
VECTOR	pDONR223
NCBI SYMBOL	CAPS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaggaggatgatataatcattcaagaacaaatgataggctgtctcaagcaatcaagatgtgctaaaaaaaactacataaaagggtgtcgtattttgactggattgggaaatatttcaacagttggacaagggaaggaatggacttttagataaggcagattttaagcaa
gctctaaaagtgtttcacttagaagtgtctgaaaaggattttgagctgtcatggctaatctgaatgacaatggcaatggcaaggttgattatggagaattcaaacgtgtgtattatgggtgaatgaatgaatacaggaaatcatatgtcgaaaggcctttatgaaactggatttcaacaaagtgg
cagtggtcctattataacataagaaaatgttactgtgcaaagaagcatttcaagtaatttcaggccattcaacagagggaagaatcaaatcatcttctagaacattaaaagttgcctgcagcaagtcgtgaagtgtcatatgtgtgaatttgaagattactatgaaggttaagtataggaat
agtagatgatgaagactttgttaacatcttacgtactccatgggggatttag