

KIAA1486-KIAA1486 protein Gene

Reference: PTXBC131727

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1486; neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adapter 2; neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adaptor 2		

For research use only.

GENE DESCRIPTION	KIAA1486 protein
GENE ID	57624
ORIGIN SPECIES	Human
ACCESSIONS	BC131727
VECTOR	pDONR223
NCBI SYMBOL	KIAA1486
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatacctccaagatgatgattctaatcctgaggaagaccctttggacacatttccagttacattgaggataggggatgaagcctacgatggcttggtattcagaatgcgtcagatattgctcgagagaatgatcgttgagaaatgaaactaacctagcctatttgaagagagaagaat
gaaaaacgccgaagacaagaagaagcaataaagcgaatcaacagaggaactgaaagtcgaagccacagcagcggagccattaccacaaagttggacaacaaggaaagaggccaccatgggctctctccagagagcctgtcaaaagctcagggaatgggatggaacaccaggggc
acctgtgtcaccagtcgactaggaagatgctgtgagccccacctgttagcgggaaccacagttcagagcctaa