



Reference: PTXBC131728

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AARSL; COXPD8; LKENP; MT-ALARS; MTALARS; alanine--tRNA ligase, mitochondrial; alanine tRNA ligase 2, mitochondrial (putative); alanyl-tRNA synthetase 2, mitochondrial (putative); alanyl-tRNA synthetase like; alanyl-tRNA synthetase 2, mitochondrial		

For research use only.

GENE DESCRIPTION	alanyl-tRNA synthetase 2, mitochondrial (putative)
GENE ID	57505
ORIGIN SPECIES	Human
ACCESSIONS	BC131728
VECTOR	pDONR223
NCBI SYMBOL	AARS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]