

NAG-neuroblastoma-amplified protein Gene

Reference: PTXBC131735

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NAG; CMT2V; MPS-III B; MPS3B; UFHSD; alpha-N-acetylglucosaminidase; N-acetylglucosaminidase, alpha; testicular tissue protein Li 18; N-acetyl-alpha-glucosaminidase		

For research use only.

GENE DESCRIPTION

neuroblastoma-amplified protein

GENE ID

51594

ORIGIN SPECIES

Human

ACCESSIONS

BC131735

VECTOR

pDONR223

NCBI SYMBOL

NAG

LENGTH

ProteoGenix bp

SEQUENCE PRIMERS

Forward primer M13R (5'→3'): GTAAACGACGGCCAGT
Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr