

ZNF138-zinc finger protein 138 Gene

Reference: PTXBC131762

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pHZ-32; zinc finger protein 138; zinc finger protein 138 (clone pHZ-32)		
For research use only.			

GENE DESCRIPTION	zinc finger protein 138
GENE ID	7697
ORIGIN SPECIES	Human
ACCESSIONS	BC131762
VECTOR	pDONR223
NCBI SYMBOL	ZNF138
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcataaattttcaaatcaaatagacacaagataagacatactgaaaaataacatttcagatgtaagaatgtgacaaatcatttgcattgttcacgcctaactcaacataaaaaattcactagagagaatttcaaaatgtgaagagtgtgaaaaaccttaactgtccacaaacctt
tctaaacctagaaaaattcactggagaaaaacctacaaatgtgaagtgtggaagcctttcaccaatcctcaatccttactaaacataagataattcgtactggagaaaaacctataaatgtcacactgtggcaagcctttaacagtcctcacaccttactagacataagataattc
atactgaagagaaacctacaaatgtgaacaatgtggcaaggtctttaagcagtcctcccaaccttactaaacatcagataatttatactggagaggaaccatacaaatgtgaggaatgtggcaagcctttaacctatcttaa