



Reference: PTXBC131763

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	matrilin-4		
For research use only.			

GENE DESCRIPTION	matrilin 4
GENE ID	8785
ORIGIN SPECIES	Human
ACCESSIONS	BC131763
VECTOR	pDONR223
NCBI SYMBOL	MATN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]