

MED1-mediator complex subunit 1 Gene

Reference: PTXBC131783

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRSP1; CRSP200; DRIP205; DRIP230; PBP; PPARBP; PPARGBP; RB18A; TRAP220; TRIP2; mediator of RNA polymerase II transcription subunit 1; ARC205; PPAR-binding protein; PPARG binding protein; TR-interacting protein 2; TRIP-2; activator-recruited cofactor 205 k		

For research use only.

GENE DESCRIPTION	mediator complex subunit 1
GENE ID	5469
ORIGIN SPECIES	Human
ACCESSIONS	BC131783
VECTOR	pDONR223
NCBI SYMBOL	MED1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr