



Reference: PTXBC131824

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein tweety homolog 3; hTTY3; tweety homolog 3; tweety family member 3		
For research use only.			

GENE DESCRIPTION	tweety homolog 3 (Drosophila)
GENE ID	80727
ORIGIN SPECIES	Human
ACCESSIONS	BC131824
VECTOR	pDONR223
NCBI SYMBOL	TTYH3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]