



Reference: PTXBC132684

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD126; IL-6R-1; IL-6RA; IL6Q; IL6RA; IL6RQ; gp80; interleukin-6 receptor subunit alpha; CD126 antigen; IL-6 receptor subunit alpha; IL-6R 1; membrane glycoprotein 80; interleukin 6 receptor		

For research use only.

GENE DESCRIPTION	interleukin 6 receptor
GENE ID	3570
ORIGIN SPECIES	Human
ACCESSIONS	BC132684
VECTOR	pDONR223
NCBI SYMBOL	IL6R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr