

C3orf35-chromosome 3 open reading frame 35 Gene

Reference: PTXBC132691

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APRG1; AP20 region protein 1; AP20 region protein1; APRG1 tumor suppressor candidate; chromosome 3 open reading frame 35		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 35
GENE ID	339883
ORIGIN SPECIES	Human
ACCESSIONS	BC132691
VECTOR	pDONR223
NCBI SYMBOL	C3orf35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaacaatggcaacaagaaaaaggtgtaaacctccagaacaggccagaatttgaataatgtgataaagaggttattgtgcgccgaactttcacacaagaattggaggagatttaacacatggaataataaacagaggaaggcgggctaagcagagcagatgggcctgcaggc
agtgctcagcatttcaacatcttccctggaccttggacccaagggtgattgcatgctggtgtctttagctaaaaataaagttaatgtcttgggcaatcataaatatggctgctcagtttctggggtgattggtggtctaaagtttcaagaacctactatgtgaaaggcatttga