

LOC440905-hypothetical protein LOC440905 Gene

Reference: PTXBC132694

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC440905		

For research use only.

GENE DESCRIPTION	hypothetical protein LOC440905
GENE ID	440905
ORIGIN SPECIES	Human
ACCESSIONS	BC132694
VECTOR	pDONR223
NCBI SYMBOL	LOC440905
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagaatttcttaattgaactctggttaatccaaaagtgtcaatcactatgtgtgggaagtttctgtacaaggaaaaaaaaaaacctaaagtcagtgtagtctaccactgtacatctgtaacctcaatccctgaaccgggcaaaatgggttcaggtcttgcaaccttga
aattccaattccatttgagagacclttgacgaggccatagctgatttcaccaccagcaactcagaaccagctactggaatgccatcagccagcagggccctgccatctatgacttctatctgtgctcactggaaggaaaccagctaccgaaggaatataccctcatcaaccaattt
acaaatggagaatagaagttaa