

FLJ33360-FLJ33360 protein Gene

Reference: PTXBC132707

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLJ33360 protein		
For research use only.			

GENE DESCRIPTION	FLJ33360 protein
GENE ID	401172
ORIGIN SPECIES	Human
ACCESSIONS	BC132707
VECTOR	pDONR223
NCBI SYMBOL	FLJ33360
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacggcacctgcatggagagaagacggcaatggatcgtgatgatctaccttgcctttaaaggaagaaatcctccatccaacacgtagatggccctggaggatgtatgccaaagtgcagaagaacgagaaggacaaatgctacatgatctcctctatggccctgcaac
agccacatcctcttgactgaagacattggccttgccttctacaggttccctctctctgggggttctctgctctctctcaagacacatggaacatgtacacgtgcccaagatggatcgccgctcagccgagcctcccagtgctgctgcccaggctggaggacccgtctctagggaca
ttgcaatccccctgctctctgagctga