



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TMEM218-transmembrane protein 218 Gene

Reference: PTXBC132716

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

transmembrane protein 218

For research use only.

GENE DESCRIPTION	transmembrane protein 218
GENE ID	219854
ORIGIN SPECIES	Human
ACCESSIONS	BC132716
VECTOR	pDONR223
NCBI SYMBOL	TMEM218
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctggtccactgtgcctcggagtcgggtcgggcgtgtcatcttagccctgctcgggtggcagtgctgctgctgtgctgctgctgcacagaccctccggggcggcagggttctctgcatctttttaattctcgtgctgtgaatcaaacatcagtctgtgctttcccgcgagctggfagaattccacccccagaagtgaaattgaattgtgatgaacttttcatgtgccctatgacctctgctgcttccatgaagccalccttctggaggccctctcttgttttaalccattatgtctcggaccgaictatgccaaacacactgcactcctaactga