



Reference: PTXBC132727

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LTrpC4; PFHB1B; TRPM4B; hTRPM4; transient receptor potential cation channel subfamily M member 4; calcium-activated non-selective cation channel 1; long transient receptor potential channel 4; melastatin-4		

For research use only.

GENE DESCRIPTION	transient receptor potential cation channel, subfamily M, member 4
GENE ID	54795
ORIGIN SPECIES	Human
ACCESSIONS	BC132727
VECTOR	pDONR223
NCBI SYMBOL	TRPM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]