

SLC25A30-solute carrier family 25, member 30 Gene

Reference: PTXBC132739

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KMCP1; kidney mitochondrial carrier protein 1; solute carrier family 25 member 30		
For research use only.			

GENE DESCRIPTION	solute carrier family 25, member 30
GENE ID	253512
ORIGIN SPECIES	Human
ACCESSIONS	BC132739
VECTOR	pDONR223
NCBI SYMBOL	SLC25A30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagccctcaactggaagccgtttgtgtacgggggctggcctccatcactgctgagtcggttacattccaattgatttaaccaagacagggctccagattcaaggccagacgaatgatgcaaaatttaaggaaattagataccagggaattgtgcacgcattagttaggtaggcaga
gaagaagggtgaaagcactctactcggggattgccccgcgatgttacgccagggcatcctatggcaccatcaagataggcacttacagagcttgaagcgactattcattgaacgccagaagatgaaactctaccgataaatgtgatgtggaattctgtctggagtcatactcaacca
ttgtaataccaactgatgtttgaaattcggatgcaagcgcaaaagcaacaccattcaaggaggaatgataggcaacttcattgaacattaccagcaagaggggacaagagagctgtggaagggtgtgtcccttactcgcagagggctgctattgtttgtgtgagctccggtctatg
acatcaccaagaagcatcttattctcagggcctgatgggagacactgtgtatacccaactcctcctcaagcttcacgtgtgtgtgtgagggggcctgtggtgtgtgagacacgtatgatgaatcagagagtgcttcagatggcagatgttctggctacacaggaacc
ctggattgctgtacagacatggaagaatgaagggtttttgtctctataaaggcttttggccaaattggttgagacttggtccttggaatatcattttctgtgacatacagcagttgaagaattgatttga