

FLJ45537-FLJ45537 protein Gene

Reference: PTXBC132741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C9orf170; chromosome 9 open reading frame 170		

For research use only.

GENE DESCRIPTION	FLJ45537 protein
GENE ID	401535
ORIGIN SPECIES	Human
ACCESSIONS	BC132741
VECTOR	pDONR223
NCBI SYMBOL	FLJ45537
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggagcaggcggcctcggcgtcagcagagcctccctgcatctcctctgggggtgtggggccgagtgaggagaactggcggacagagaaaagggcctcgtcgcgaagcccgccgggggagcctggcctcctgctcagtcggtgcaaatgaaaaagagacgttctgttctccgaaaacctgactaatacagtagaggatccagattgacaactcaggagaaagtctgattaggagtcggaagccagattgaaaaacctctgatagatgtcactcgggaagaccatgaaaacagccaaataatcaagagaagatgcaaatgaaactgtgaacagaccaaagatga