

LOC389791-hypothetical LOC389791 Gene

Reference: PTXBC132745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	hypothetical LOC389791
GENE ID	389791
ORIGIN SPECIES	Human
ACCESSIONS	BC132745
VECTOR	pDONR223
NCBI SYMBOL	LOC389791
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgtgaagacgagagcgaagcgaggacgacgctcggggaagcgggtgacccgaggagactgctggagcgaactggcgttccgggggctcctccgggaaggggaatcgggatgtggctcgaaggaccgaaggcccacctcgaccgccccgagtggt
ttgagttgtagtgctgtctcggatgcagggcgagcaccgagcgggtaaccagtcctcctgagagctgcaggtcccagccccgttttacagatagggaaactgaggcagctgctgggacatcgccataggggtttgctccaccactcctcagccactttacagaccaggag
ttgaggactgctcgaacttagaggttttaggggtgtgagggtgacccccatgtctgtgctctactatccctatccaggccttcccagcggcactctcagctcctgctgtagctcctgctgtagctcctgctcctccagcttcctcctctgtttacattaattgctcctgatctt
ccacacttga