



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C13orf39-chromosome 13 open reading frame 39 Gene

Reference: PTXBC132748

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf39; protein-lysine methyltransferase METTL21C; methyltransferase-like protein 21C; methyltransferase like 21C		

For research use only.

GENE DESCRIPTION	chromosome 13 open reading frame 39
GENE ID	196541
ORIGIN SPECIES	Human
ACCESSIONS	BC132748
VECTOR	pDONR223
NCBI SYMBOL	C13orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgtgtgtctgagctccgcgcagcagcctggcgccgcgggggaagcactcagctcccgcgggtgctgctgttagagggctgagaagaagggggctccgcagaaagacagcaccgggggagtcctagaagaatccaaacagatagaaccactctctcatagcctccagaattgtt
cctcagattacgccgactacactcaggagcattatcggtttgcaggaaaaggagattgtatccaggaaatccatagagagttacggagcgggtgtgtgtggccgggggctatggcttctgtcctaacttggaggaaactgccgagggaattgaattccaaagatcgaataaactgtgaattggt
ggcggaaccaggccgtgttccattgtggccgagatttttagagctcgaagctacagacagagattgctgatgctcctggaaaccttcaatacaacttttaaaaaaacacactacaatgtacacagacatctgctcgaagtgaagaacctggatgggggaagacctggcaaaaaacttcccaa
gtcagcttttatactgattacgtcctagctcagatgtgtgtaccactcactacttccggaacagctccaccacatggtgtactcttccagcggggagctgctgtctgtgtggcaaacaaattcaggttcagccggactatgaattttatagaattcaagcaagttttgacacaacact
gttgcctcaataccagagcttccagctcaaaacttttlaaggggatcctacaaatgagactaa

More info: Please contact us for any question: contact@proteogenix.fr