

C12orf69-chromosome 12 open reading frame 69 Gene

Reference: PTXBC132752

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf69; single-pass membrane and coiled-coil domain-containing protein 3; single-pass membrane protein with coiled-coil domains 3		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 69
GENE ID	440087
ORIGIN SPECIES	Human
ACCESSIONS	BC132752
VECTOR	pDONR223
NCBI SYMBOL	C12orf69
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaagtgactcttaccagagaaacccaaaagcgaggagaagtaaatcgttcaccagcagctcttgattgcttatctgacagcttcgatgcacaaatagctgactgaggttctaataatgcactgggtgcaggctggcctcattgagatgaaaagagatgggacca
tcaagaaaactgtgacctcatcatcgaagccattatgaaatcctaaaggaattgcagaaggttgatgaagcactaaaagataagctagagccaacctctatagaaaactcaggatattaaggaaaaggaaacagacaaaatgcaatagtgcaaaagggttattcggtcatctgggag
aagctacatctgcagccagtcagtcgctgttaaaactgtgggctcaaatgtcacaactggcacaataaacaagttggtcactgtgttagctcaaattggtgcttctccttggttagtattggagttgctgttcttgccctggcatagatgattgtccgtgccatcctgggagcagtggaataaa
cacagctcaagcagccatcaaaagttagagaagcatctgggtgaggtcaaatcagcctcagaaaataataatcatgccattactgaggtcatcaatacagtgaaacaccaaataaatga