

LOC123688-similar to RIKEN cDNA C630028N24 gene Gene

Reference: PTXBC132753

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to RIKEN cDNA C630028N24 gene
GENE ID	123688
ORIGIN SPECIES	Human
ACCESSIONS	BC132753
VECTOR	pDONR223
NCBI SYMBOL	LOC123688
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaagtggaaactatcagcagtcagaggtcttagcaaacctttcagtgaggaaacgcctctgcttagtgaggcagtcgtttgggtgaaagtccagggtccgacctctagctatgatgacaaaacttcatgtctacgtttcaaaaaccaaagatggcccaactgaatat
gtcctcaaaaataagcaacaccaaaggctagcaaaaatccagacctgattgaagtgcagaatcacatcatatgtttctgaaagccgtggatttccaacagcctctgtgtgtcacactaaaggagacaacacagcttctctgtgtgtatagtggtctgaaatcaaaagctacttgggtgag
gctgctgacttacctccaggaagacctatcgtgagcttccgtcagcccccagctattgtatgaattggaaaactagctgccaaattggataagacactgcagagattccatcacccaaagttaagttagtcttcacgggagaactcatctggaatctgaaaaatgttctctctggagaa
atacctgtatgcctggccagaatcgaaaccgagagattgtgagcatgtcattcatctgttcaaggaggaagttaaccaaattaagtcatttgcagaatgtgagtatttcccaattaa