

LOC401296-hypothetical protein LOC401296 Gene

Reference: PTXBC132760

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC401296
GENE ID	401296
ORIGIN SPECIES	Human
ACCESSIONS	BC132760
VECTOR	pDONR223
NCBI SYMBOL	LOC401296
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactggggaggagggaagcactctcactctgctggagcctctgccaggaaaccccgattggccagccgctccctggggaggcccatcccaagctgtcccgaccccggtgtctgtcagcactgatcagggtctccagacactgtctgtgggtgcatggggacgagcgggagggtgagtggtggcgacacagctgtctgtagacacacagaggcatgcacgcatgtgtctgtgtcactgcatgtgttcacacgactgggcacccggcggagggtttcacacagcgggtgtacaccgacccctctgcagatgtgcacacccagcaaggctacgtccacaggcgtgtatgtttctgcaccagggaagtgggtcacgccccacaccacatcgtatactctgtccctgcggcacgcacctgtctgggtgcaccaggcagcttaagcccatgtcacctatgtgtacacatgtgtgcacatactactcaggcacaccagtggtgtgtacacggctgtccaggcacactgggcagttcaagccacgcctgccagactctggaatcgggcactttctga