

FLJ45256-hypothetical LOC400511 Gene

Reference: PTXBC132761

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC400511
GENE ID	400511
ORIGIN SPECIES	Human
ACCESSIONS	BC132761
VECTOR	pDONR223
NCBI SYMBOL	FLJ45256
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttagtagacgggtttgccatgttgccagctggtctcgaactcctgagctcaagtgaagggccggcctcgacctccaaagtctgggattacaggcatgagccaccgcgccttgccctttgctgatttttaacccgcatcctccactgtgtaagctgtaactgtgaata
caacagctttttgagttctgcgattctttgagcaaatgctacactgagcgtgcttgggactcccgacacactctccactccctctctgtttctcactctctctggagtaaacgaaagaggccctgcttgcctcttgcaggcagttctgaacacactgagatcattgtgatctcatac
ctcttctaaagcaatcatcaagaattctttcaataaaccagtgctcttcattggaagtgatgctgatgaaatgactttccagcttattacagtttcagaaagctcctgtttaagagtttgagcatgtatgtgtgtgcatatacaccagcacacacatgcatgccacagctgcatgccttc
accagaaccaggacgaggagtattctattgccagaattga